



# Anti-*P. aeruginosa* Tail Sheath Polyclonal antibody (DPAB-DC4676)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Target	<i>P. aeruginosa</i> Tail Sheath
Immunogen	A synthetic peptide corresponding to internal region of <i>Pseudomonas aeruginosa</i> Tail Sheath.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	<i>P. aeruginosa</i>
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Introduction	<i>Pseudomonas aeruginosa</i> is a common bacterium that can cause disease in animals, including humans. It is citrate, catalase, and oxidase positive. It is found in soil, water, skin flora, and most man-made environments throughout the world. It thrives not on
--------------	--

**Keywords**

Tail Sheath; Pseudomonas aeruginosa Tail Sheath

---