



Anti-*P. aureofaciens* Polyclonal antibody (DPAB-DC4678)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	<i>P. aureofaciens</i>
Immunogen	<i>Pseudomonas aureofaciens</i>
Source/Host	Rat
Species Reactivity	<i>P. aureofaciens</i>
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Lyophilized
Size	100 µl
Preservative	None
Storage	Store at 4°C on dry atmosphere. After reconstitution with 50 uL of distilled water and 50 uL of glycerol, store at -20°C for 2 years. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	<i>Pseudomonas aureofaciens</i> is a yellowish, aerobic, Gram-negative, motile, polar flagellated, rod bacterium isolated from clay near the river Maas. Based on 16S rRNA analysis, <i>P. aureofaciens</i> has been placed in the <i>P. chlororaphis</i> group.
Keywords	<i>P. aureofaciens</i> ; <i>Pseudomonas aureofaciens</i>