



Anti-*P. marginalis* Polyclonal antibody (DPAB-DC4679)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|--------------------|---|
| Target | <i>P. marginalis</i> |
| Immunogen | <i>Pseudomonas marginalis</i> |
| Source/Host | Rat |
| Species Reactivity | <i>P. marginalis</i> |
| Conjugate | Unconjugated |
| Applications | WB, ELISA, |
| Format | Lyophilized |
| Size | 100 µl |
| Preservative | None |
| Storage | Store at 4°C on dry atmosphere. After reconstitution with 50 µL of distilled water and 50 µL of glycerol, store at -20°C for 2 years. Aliquot to avoid repeated freezing and thawing. |

BACKGROUND

| | |
|--------------|---|
| Introduction | <i>Pseudomonas marginalis</i> is a soil bacterium that can cause soft rots of plant tissues. It infects poinsettia, lettuce, crucifers (canola, mustard). |
| Keywords | <i>P. marginalis</i> ; <i>Pseudomonas marginalis</i> |