



Anti-Ractopamine polyclonal antibody (DPAB-DC4684)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Ractopamine conjugated with BSA.
Source/Host	Sheep
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	EIA,
Format	Liquid
Size	100 µl
Buffer	In 20 mM phosphate buffer, 150 mM NaCl, pH 7.2 (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Ractopamine is a drug used as a feed additive to promote leanness in animals raised for their meat. Pharmacologically, it is a beta-adrenergic agonist. It is the active ingredient in products known as Paylean for swine and Optaflexx for cattle, developed
Keywords	Ractopamine