



Anti-RAF1 polyclonal antibody (DPAB-DC4686)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to human RAF1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Purification	Protein A purification
Conjugate	Unconjugated
Applications	WB (Cell lysate),
Format	Liquid
Concentration	1 mg/mL
Size	100 µl
Buffer	In HEPES, 150 mM NaCl (50% glycerol, 0.01% BSA, 0.03% sodium azide)
Preservative	0.03% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	RAF1
-----------------	------

GENE INFORMATION

Entrez Gene ID [5894](#)

UniProt ID [P04049](#)
