



# Anti-RanGAP polyclonal antibody (DPAB-DC4687)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	Reacts with human & Xenopus RanGAP.
<b>Immunogen</b>	A mixture of synthetic peptide corresponding to amino acids and The sequence is CHWSDMFTGRLRPEI, CPSPEKLVRMGPRRSA
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Clawed frog, Frog, Human
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, IF,
<b>Format</b>	Liquid
<b>Size</b>	100 µl
<b>Buffer</b>	In PBS (0.1% sodium azide)
<b>Preservative</b>	0.1% Sodium Azide
<b>Storage</b>	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

<b>Keywords</b>	RanGAP
-----------------	--------