



# Anti-RBC (native) polyclonal antibody (DPAB-DC4689)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Native purified human Red Blood Cells (RBC).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	Hemagglutinin Inhibition,
Format	Lyophilized
Size	2 ml
Buffer	In 10 mM sodium phosphate buffer, 140 mM NaCl, pH 7.4 (0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C on dry atmosphere. After reconstitution with 0.1 mL of deionized water, store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Keywords	RBC
----------	-----