



# Anti-RCC1 polyclonal antibody (DPAB-DC4690)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                           |   |
|---------------------------|---|
| <b>Specificity</b>        | Reacts with human & Xenopus RCC1.   |
| <b>Immunogen</b>          | Recombinant His fusion protein corresponding to Xenopus RCC1.                                       |
| <b>Isotype</b>            | IgG   |
| <b>Source/Host</b>        | Rabbit  |
| <b>Species Reactivity</b> | Clawed frog, Frog, Human  |
| <b>Conjugate</b>          | Unconjugated  |
| <b>Applications</b>       | WB (Cell lysate), IF, IP,   |
| <b>Format</b>             | Liquid  |
| <b>Size</b>               | 200 µl  |
| <b>Buffer</b>             | In PBS (0.1% sodium azide)  |
| <b>Preservative</b>       | 0.1% Sodium Azide   |
| <b>Storage</b>            | Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing. |

## BACKGROUND

|                 |      |
|-----------------|------|
| <b>Keywords</b> | RCC1 |
|-----------------|------|

## GENE INFORMATION

Entrez Gene ID [1104](#)

---

UniProt ID [P18754](#)

---