



Anti-RSV Matrix protein 2 Polyclonal antibody (DPAB-DC4696)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	RSV Matrix protein 2
Immunogen	A synthetic peptide corresponding to amino acids 81-95 of human RSV M2 protein. The sequence is LESYIGSINNITKQSA
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	RSV
Conjugate	Unconjugated
Applications	WB (Cell lysate),
Format	Liquid
Size	100 µg
Buffer	In PBS (0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Human respiratory syncytial virus (RSV) is a virus that causes respiratory tract infections. It is a major cause of lower respiratory tract infections and hospital visits during infancy and
--------------	---

childhood. A prophylactic medication (not a vaccine) exists for p

Keywords

RSV M2 protein; Respiratory syncytial virus
