



## Anti-RSV RSS-2 Polyclonal antibody (DPAB-DC4698)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Specificity</b>	human respiratory syncytial virus (RSV) Fusion protein of RSS-2 (subtype A) strain.
<b>Target</b>	RSV RSS-2
<b>Immunogen</b>	A synthetic peptide corresponding to amino acids 201-211 of human respiratory syncytial virus. The sequence is KQLLPIVNKQSC
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	RSV
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB,
<b>Format</b>	Liquid
<b>Size</b>	100 µg
<b>Buffer</b>	In PBS (0.02% sodium azide)
<b>Preservative</b>	0.02% Sodium Azide
<b>Storage</b>	Store at 4°C for short term. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

### BACKGROUND

**Introduction**

Human respiratory syncytial virus (RSV) is a virus that causes respiratory tract infections. It is a major cause of lower respiratory tract infections and hospital visits during infancy and childhood. A prophylactic medication (not a vaccine) exists for p

---

**Keywords**

RSV RSS-2; Respiratory syncytial virus RSS-2

---