



Anti-*S. cerevisiae* RFA1/RFA2/RFA3 Polyclonal antibody (DPAB-DC4700)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	<i>S. cerevisiae</i> RFA1/RFA2/RFA3
Immunogen	Recombinant protein corresponding to <i>Saccharomyces cerevisiae</i> RFA1/RFA2/RFA3.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	<i>S. cerevisiae</i>
Conjugate	Unconjugated
Applications	WB, IP,
Format	Lyophilized
Size	100 µl
Buffer	Lyophilized from serum
Preservative	None
Storage	Store at -20°C on dry atmosphere. After reconstitution with 100 uL deionized water, store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	RFA1/RFA2/RFA3
----------	----------------