



Anti-Ritodrine polyclonal antibody (DPAB-DC4704)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Ritodrine conjugated with BSA.
Source/Host	Sheep
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	WB, EIA,
Format	Liquid
Size	100 µl
Buffer	In 20 mM phosphate buffer, 150 mM NaCl, pH 7.2 (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Ritodrine (discontinued preparation: Yutopar) is a tocolytic drug, used to stop premature labor. This drug has been removed from the US market, according to FDA Orange Book. It was available in oral tablets or as an injection and was typically used as the
Keywords	Ritodrine