



Anti-Ronidazole polyclonal antibody (DPAB-DC4705)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Ronidazole conjugated with BTG.
Source/Host	Sheep
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	EIA,
Format	Liquid
Size	100 µl
Buffer	In 20 mM phosphate buffer, 150 mM NaCl, pH 7.2 (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Ronidazole is an antiprotozoal agent used in veterinary medicine for the treatment of histomoniasis and swine dysentery. It may also have use for the treatment of Tritrichomonas foetus infection in cats.
Keywords	Ronidazole