



Anti-RORB (aa 111-459) polyclonal antibody (DPAB-DC4706)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant LBD fusion protein corresponding to amino acids 111-459 of RORB.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse
Conjugate	Unconjugated
Applications	WB (Cell lysate), WB (Recombinant protein), IF, GSA,
Format	Liquid
Size	100 µl
Buffer	In PBS (0.01% thimerosal)
Preservative	None
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	RORB
----------	------

GENE INFORMATION

Entrez Gene ID [6096](#)

UniProt ID [Q58EY0](#)
