



# Anti-Rotavirus Polyclonal antibody (DPAB-DC4707)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	ICPs and late structural (virion) antigens. Cross reactivity is >90% with human rotaviruses (reported). Uninfected cell reactivity is negative against HEP-2 cells and WI-38 cells by indirect immunofluorescence.
<b>Target</b>	Rotavirus
<b>Immunogen</b>	Bovine (Nebraska Calf Diarrhea Virus).
<b>Source/Host</b>	Goat
<b>Species Reactivity</b>	Rotavirus
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IF,
<b>Format</b>	Liquid
<b>Size</b>	1 ml
<b>Buffer</b>	In 10 mM PBS, pH 7.2 (0.1% sodium azide)
<b>Preservative</b>	0.1% Sodium Azide
<b>Storage</b>	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

**Keywords** Rotavirus NCDV