



Anti-drosophila Rtf1 Polyclonal antibody (DPAB-DC4710)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------------------|---|
| Target | drosophila Rtf1 |
| Immunogen | Recombinant protein corresponding to full length drosophila Rtf1. |
| Isotype | IgG |
| Source/Host | Rabbit |
| Species Reactivity | drosophila |
| Conjugate | Unconjugated |
| Applications | WB (Cell lysate), IF, |
| Format | Liquid |
| Size | 50 µg |
| Buffer | In PBS (0.05% sodium azide, 0.05% proclin 300) |
| Preservative | 0.05% Sodium Azide |
| Storage | Store at -20°C. For long term storage store at -80°C. Aliquot to avoid repeated freezing and thawing. |

BACKGROUND

| | |
|---------------------|---|
| Introduction | Drosophila is a genus of small flies, belonging to the family Drosophilidae, whose members are often called "fruit flies" or (less frequently) pomace flies, vinegar flies, or wine flies, a reference to the characteristic of many species to linger around o |
|---------------------|---|

