



Anti-Rubella virus Virions Polyclonal antibody (DPAB-DC4711)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Purified virions. Negative against VERO cells by direct immunofluorescence.
Target	Rubella virus Virions
Immunogen	Rubella Virus, strain HPV77.
Source/Host	Goat
Species Reactivity	Rubella virus
Conjugate	Unconjugated
Applications	IF, Hemagglutinin Inhibition,
Format	Liquid
Size	1 ml
Buffer	In 10 mM PBS, pH 7.2 (0.1% sodium azide)
Preservative	0.1% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Rubella virus is the pathogenic agent of the disease Rubella, and is the cause of congenital rubella syndrome when infection occurs during the first weeks of pregnancy. Humans are the only known host of this virus.
--------------	--

Keywords

Virions; Rubella Virus Virions
