



# Anti-RUNX1T1 polyclonal antibody (DPAB-DC4712)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Specificity	This antibody recognizes the ETO (RUNX1T1) part of the fusion protein.
Immunogen	A synthetic peptide (conjugated with KLH) corresponding to RUNX1T1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Cell lysate), WB (Recombinant protein), ELISA,
Format	Liquid
Size	50 µg
Buffer	In PBS (0.05% sodium azide, 0.05% proclin 300)
Preservative	0.05% Sodium Azide
Storage	Store at -20°C. For long term storage store at -80°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Keywords	RUNX1T1
----------	---------

# GENE INFORMATION

Entrez Gene ID [862](#)

---

UniProt ID [Q06455](#)

---