



Anti-RVSV Nucleocapsid Polyclonal antibody (DPAB-DC4713)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	RVSV Nucleocapsid
Immunogen	A synthetic peptide corresponding to 17 amino acids at internal region of the RVF virus nucleocapsid.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	RVSV
Purification	Peptide affinity purification
Conjugate	Unconjugated
Applications	ELISA,
Format	Liquid
Concentration	1 mg/mL
Size	100 µg
Buffer	In PBS (0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C for three months. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction

Rift Valley fever (RVF) is a viral disease that can cause mild to severe symptoms. The mild symptoms may include: fever, muscle pains, and headaches which often last for up to a week. The severe symptoms may include: loss of the ability to see beginning t

Keywords

RVFV Nucleocapsid
