



Anti-RVFP Polyclonal antibody (DPAB-DC4714)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	RVFV
Immunogen	A synthetic peptide corresponding to internal region 20 amino acids of Rift Valley Fever Virus.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	RVFV
Conjugate	Unconjugated
Applications	ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS (0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C for three months. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Rift Valley fever (RVF) is a viral disease that can cause mild to severe symptoms. The mild symptoms may include: fever, muscle pains, and headaches which often last for up to a week. The severe symptoms may include: loss of the ability to see beginning t
Keywords	RVFV