



Rabbit anti-Salmonella Salmonella polyclonal antibody (DPAB-DC4717)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Polyvalent for Salmonella "O" & "H" antigens. Immunocaptures Salmonellae. Antiserum is not absorbed and does react with related Enterobacteriaceae.
Immunogen	A mixture of Salmonella enteritidis, Salmonella typhimurium, and Salmonella heidelberg.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Salmonella
Conjugate	Unconjugated
Applications	IF, ELISA,
Format	Liquid
Size	100 µl
Buffer	In 10 mM PBS, pH 7.2 (0.1% sodium azide)
Preservative	0.1% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Salmonella /s?lm??n?l?/ is a genus of rod-shaped, Gram-negative bacteria. There are only two species of Salmonella, Salmonella bongori and Salmonella enterica, of which there are
---------------------	--

around six subspecies and innumerable serovars. The genus belongs to the sa

Keywords

Salmonella
