



Anti-*S. mediterranea* Dvl-b Polyclonal antibody (DPAB-DC4723)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	<i>S. mediterranea</i> Dvl-b
Immunogen	A synthetic peptide corresponding to C-terminus of Freshwater Flatworms <i>Schmidtea mediterranea</i> Dvl-b.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	<i>S. mediterranea</i>
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	<i>Schmidtea mediterranea</i> is a freshwater planarian that lives in southern Europe and Tunisia. It is a model for regeneration and development of tissues such as the brain and germline.
---------------------	--

Keywords

Dvl-b; Schmidtea mediterranea Dvl-b
