



Anti-IAV H1N1 Hemagglutinin Polyclonal antibody (DPAB-DC4725)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	The peptide sequence is unique from the peptide sequence for product DPAB-DC4724 and DPAB-DC4726. This antibody is a cognate pair with antibody DPAB-DC4799.
Target	IAV H1N1 Hemagglutinin
Immunogen	A synthetic peptide corresponding to Seasonal H1N1 Hemagglutinin.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	IAV
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS (0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C for three months. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction

Influenza A (H1N1) virus is the subtype of influenza A virus that was the most common cause of human influenza (flu) in 2009. It is an Orthomyxovirus that contains the glycoproteins Haemagglutinin and Neuraminidase. For this reason, they are described as

Keywords

Hemagglutinin; Seasonal H1N1 Hemagglutinin
