



Anti-IAV Nonstructural Protein 1 Polyclonal antibody (DPAB-DC4729)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	IAV Nonstructural Protein 1
Immunogen	A synthetic peptide corresponding to Seasonal Influenza NS1 protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	IAV
Conjugate	Unconjugated
Applications	ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS (0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C for three months. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Influenza, commonly known as "the flu", is an infectious disease of birds and mammals caused by RNA viruses of the family Orthomyxoviridae, the influenza viruses. The most common symptoms are chills, fever, runny nose, sore throat, muscle pains, headache
---------------------	--

Keywords

NS1 protein; Seasonal Influenza NS1 protein
