



# Anti-Secobarbital polyclonal antibody (DPAB-DC4730)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	Phenobarbital <23.1%, Pentobarbital <23.1%, Butalbital <23.1%, Amobarbital <23.1%, Barbital <23.1%
<b>Immunogen</b>	Secobarbital conjugated with BSA.
<b>Source/Host</b>	Sheep
<b>Species Reactivity</b>	N/A
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	EIA,
<b>Format</b>	Liquid
<b>Size</b>	1 ml
<b>Buffer</b>	In 20 mM phosphate buffer, 150 mM NaCl, pH 7.2 (0.09% sodium azide)
<b>Preservative</b>	0.09% Sodium Azide
<b>Storage</b>	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

<b>Introduction</b>	Secobarbital sodium (marketed by Eli Lilly and Company, and subsequently by other companies as described below, under the brand name Seconal) is a barbiturate derivative drug that was patented in 1934 in the US. It possesses anaesthetic, anticonvulsant, a
<b>Keywords</b>	Secobarbital