



Anti-5-HT polyclonal antibody (DPAB-DC4731)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Serotonin conjugated with BSA.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Rat
Conjugate	Unconjugated
Applications	IHC-P, IHC-Fr,
Format	Lyophilized
Size	50 µl
Buffer	Lyophilized from serum
Preservative	None
Storage	Store at 4°C on dry atmosphere. After reconstitution with 50-100 uL distilled water, store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Serotonin or 5-hydroxytryptamine (5-HT) is a monoamine neurotransmitter. Biochemically derived from tryptophan, serotonin is primarily found in the gastrointestinal tract (GI tract), platelets, and the central nervous system (CNS) of animals, including hu
Keywords	5-HT; Serotonin (5-HT)