



Anti-5-HT polyclonal antibody (DPAB-DC4732)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Hapten conjugated Serotonin (5-Hydroxy-tryptamine).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	IHC-P, IHC-Fr,
Format	Liquid
Size	1 ml
Buffer	In phosphate buffer (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at 4°C. Do not freeze.

BACKGROUND

Introduction	Serotonin or 5-hydroxytryptamine (5-HT) is a monoamine neurotransmitter. Biochemically derived from tryptophan, serotonin is primarily found in the gastrointestinal tract (GI tract), platelets, and the central nervous system (CNS) of animals, including hu
Keywords	5-HT; Serotonin