



Anti-SFTSV HB29 (C-terminal) polyclonal antibody (DPAB-DC4737)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	SFTSV HB29 Membrane Glycoprotein antibody is specific to the SFTSV HB29.
Immunogen	A synthetic peptide corresponding to 15 amino acids at C-terminus of human SFTSV HB29.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Virus
Purification	Peptide affinity purification
Conjugate	Unconjugated
Applications	ELISA,
Format	Liquid
Concentration	1 mg/mL
Size	100 µg
Buffer	In PBS (0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C for three months. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction

SFTSV is a phlebovirus in the family Bunyaviridae. It appears to be more closely related to the Uukuniemi virus serogroup than to the Sandfly fever group. It is a member of the Bhanja virus serocomplex. The clinical condition it causes is known as severe

Keywords

SFTSV HB29
