



Anti-Shaker (N-terminal) polyclonal antibody (DPAB-DC4738)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody recognizes ~73 KDa of Shaker protein.
Immunogen	A synthetic peptide corresponding to N-terminus of Shaker.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Fruit fly, Human
Conjugate	Unconjugated
Applications	WB, IHC, IP, ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS, pH 7.2 (10% Proclin300)
Preservative	None
Storage	Store at 4°C. For long term storage store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	Shaker
----------	--------