



Anti-Sir2.1 (full length) polyclonal antibody (DPAB-DC4742)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to full length sir-2.1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Nematoda
Conjugate	Unconjugated
Applications	WB (Recombinant protein), IHC,
Format	Liquid
Size	100 µl
Buffer	In antiserum (0.05% sodium azide)
Preservative	0.05% Sodium Azide
Storage	Store at -20°C or -80°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	Sir2.1
----------	--------