



# Anti-SLC polyclonal antibody (DPAB-DC4743)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Recombinant protein corresponding to rat SLC.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Rat
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB,
<b>Format</b>	Liquid
<b>Size</b>	100 µg
<b>Buffer</b>	In PBS, pH 7.2 (50% glycerol, 1% BSA, 0.02% sodium azide)
<b>Preservative</b>	0.02% Sodium Azide
<b>Storage</b>	Store at -20°C. For long term storage store at -80°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

<b>Keywords</b>	SLC
-----------------	-----