



Anti-Somatostatin polyclonal antibody (DPAB-DC4745)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody reacts with somatostatin 1-14 and somatostatin 1-28 (GRIF, SRIF).
Immunogen	Synthetic somatostatin (UCB) conjugated with ovalbumin.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Conjugate	Unconjugated
Applications	IHC-P,
Format	Liquid
Size	100 µl
Buffer	In phosphate buffer, pH 7.3 (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at 4°C. Do not freeze.

BACKGROUND

Introduction Somatostatin (also known as growth hormone-inhibiting hormone (GHIH) or somatotropin release-inhibiting factor (SRIF)) or somatotropin release-inhibiting hormone[citation needed] is a peptide hormone that regulates the endocrine system and affects neurotr

Keywords

Somatostatin

GENE INFORMATION

Entrez Gene ID

[6750](#)

UniProt ID

[P61278](#)
