



Anti-SR-beta (aa 246-265) polyclonal antibody (DPAB-DC4750)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to amino acids 246-265 of SR beta.
Source/Host	Sheep
Species Reactivity	Dog
Conjugate	Unconjugated
Applications	WB (Cell lysate), IP,
Format	Liquid
Size	250 µg
Buffer	In PBS (0.08% sodium azide)
Preservative	0.08% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	The eukaryotic SRP receptor (termed SR) is a heterodimer of SR-alpha (70 kDa) and SR-beta (25 kDa), both of which contain a GTP-binding domain, while the prokaryotic SRP receptor comprises only the monomeric loosely membrane-associated SR-alpha homologue
Keywords	SR-beta