



Anti-SSeCKS (aa 620-1596) polyclonal antibody (DPAB-DC4751)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids 620-1596 of SSeCKS.
Source/Host	Sheep
Species Reactivity	Bovine, Rat
Conjugate	Unconjugated
Applications	WB (Tissue lysate), IHC, IP,
Format	Liquid
Size	50 µg
Buffer	In PBS (0.08% sodium azide)
Preservative	0.08% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	SSeCKS
----------	--------

GENE INFORMATION

Entrez Gene ID	9590
----------------	----------------------

UniProt ID

[Q02952](#)
