



Anti-Stanozolol polyclonal antibody (DPAB-DC4752)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	16 beta-Hydroxystanozolol 40%
Immunogen	Stanozolol conjugated with BSA.
Source/Host	Sheep
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	EIA,
Format	Liquid
Size	100 µl
Buffer	In 20 mM phosphate buffer, 150 mM NaCl, pH 7.2 (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Stanozolol, commonly sold under the name Winstrol (oral) and Winstrol Depot (intramuscular), is a synthetic anabolic steroid derived from dihydrotestosterone. It was developed by Winthrop Laboratories (Sterling Drug) in 1962, and has been approved by the
Keywords	Stanozolol