



## Anti-S. aureus Eap Polyclonal antibody (DPAB-DC4753)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Target</b>	S. aureus Eap
<b>Immunogen</b>	A synthetic peptide corresponding to C-terminus of Staphylococcus aureus Eap.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	S. aureus
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, ELISA,
<b>Format</b>	Liquid
<b>Size</b>	100 µg
<b>Buffer</b>	In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)
<b>Preservative</b>	0.01% Sodium Azide
<b>Storage</b>	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

### BACKGROUND

**Introduction** Staphylococcus aureus is a Gram-positive coccal bacterium that is a member of the Firmicutes, and is frequently found in the human respiratory tract and on the skin. It is positive for catalase and nitrate reduction. Although S. aureus is not always patho

**Keywords**

Eap; Staphylococcus aureus Eap

---