



Anti-S. aureus Staphylococcus Enterotoxin B Polyclonal antibody (DPAB-DC4755)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Staphylococcus aureus Enterotoxin B. Non-reactive with other Staph. toxins.
Target	S. aureus Staphylococcus Enterotoxin B
Immunogen	Purified Staphylococcus aureus Enterotoxin B.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	S. aureus
Conjugate	Unconjugated
Applications	ID,
Format	Liquid
Size	1 mg
Buffer	In 10 mM PBS, pH 7.2 (0.1% sodium azide)
Preservative	0.1% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction Staphylococcus aureus is a Gram-positive coccal bacterium that is a member of the Firmicutes,

and is frequently found in the human respiratory tract and on the skin. It is positive for catalase and nitrate reduction. Although *S. aureus* is not always patho

Keywords

SEB; *Staphylococcus aureus* Enterotoxin B (SEB)
