



# Anti-S. aureus ScpA Polyclonal antibody (DPAB-DC4760)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Target	S. aureus ScpA
Immunogen	A synthetic peptide corresponding to C-terminus of Staphylococcus aureus scpA.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	S. aureus
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Introduction	Staphylococcus aureus is a Gram-positive cocci bacterium that is a member of the Firmicutes, and is frequently found in the human respiratory tract and on the skin. It is positive for catalase and nitrate reduction. Although S. aureus is not always patho
--------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

**Keywords**

scpA; Staphylococcus aureus scpA

---