



# Anti-S. aureus Tail Tape Measure protein Polyclonal antibody (DPAB-DC4761)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Target</b>	S. aureus Tail Tape Measure protein
<b>Immunogen</b>	A synthetic peptide corresponding to N-terminus of Staphylococcus phage phi13 Phage Tail Tape Measure protein.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	S. aureus
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, ELISA,
<b>Format</b>	Liquid
<b>Size</b>	100 µg
<b>Buffer</b>	In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)
<b>Preservative</b>	0.01% Sodium Azide
<b>Storage</b>	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

**Keywords** Tail Tape Measure protein; Staphylococcus phage phi13 Phage Tail Tape Measure protein