



Anti-*S. maltophilia* Polyclonal antibody (DPAB-DC4762)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	<i>S. maltophilia</i>
Immunogen	<i>Stenotrophomonas maltophilia</i>
Source/Host	Rat
Species Reactivity	<i>S. maltophilia</i>
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Lyophilized
Size	100 µl
Preservative	None
Storage	Store at 4°C on dry atmosphere. After reconstitution with 50 uL of distilled water and 50 uL of glycerol, store at -20°C for 2 years. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	<i>Stenotrophomonas maltophilia</i> is an aerobic, nonfermentative, Gram-negative bacterium. It is an uncommon bacterium and human infection is difficult to treat. Initially classified as <i>Pseudomonas maltophilia</i> , <i>S. maltophilia</i> was also grouped in the genus <i>Xant</i>
Keywords	<i>S. maltophilia</i> ; <i>Stenotrophomonas maltophilia</i>