



Anti-*S. maltophilia* Polyclonal antibody (DPAB-DC4763)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	<i>S. maltophilia</i>
Immunogen	<i>Stenotrophomonas maltophilia</i> total protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	<i>S. maltophilia</i>
Purification	Ammonium sulfate precipitation
Conjugate	Unconjugated
Applications	WB,
Format	Lyophilized
Size	100 µl
Buffer	Lyophilized from PBS, pH 7.4
Preservative	None
Storage	Store at 4°C on dry atmosphere. After reconstitution with 50 uL of distilled water and 50 uL of glycerol, store at -20°C for 2 years. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction *Stenotrophomonas maltophilia* is an aerobic, nonfermentative, Gram-negative bacterium. It is

an uncommon bacterium and human infection is difficult to treat. Initially classified as *Pseudomonas maltophilia*, *S. maltophilia* was also grouped in the genus *Xant*

Keywords

S. maltophilia; *Stenotrophomonas maltophilia*
