



Anti-S. aureus Group A GAS Polyclonal antibody (DPAB-DC4765)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Type specific carbohydrate of group A Streptococcus. Does not react with other Strep groups.
Target	S. aureus Group A GAS
Immunogen	Streptococci Group A.
Source/Host	Goat
Species Reactivity	S. aureus
Conjugate	Unconjugated
Applications	ELISA, sELISA, Capture Antibody, Lateral-Flow Immunochromatographic Assay, Detection Antibody,
Format	Liquid
Size	1 mg
Buffer	In 10 mM PBS, pH 7.2 (0.1% sodium azide)
Preservative	0.1% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	S. pyogenes, also known as Group A Streptococcus (GAS), is the causative agent in a wide range of Group A streptococcal infections. These infections may be non-invasive or invasive.
--------------	---

The non-invasive infections tend to be more common and less severe. The m

Keywords

GAS; Streptococcus Group A
