



Anti-S. aureus Group B GBS Polyclonal antibody (DPAB-DC4767)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Type specific carbohydrate of group B Streptococcus. Does not react with other Strep groups.
Target	S. aureus Group B GBS
Immunogen	Streptococci Group B.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	S. aureus
Conjugate	Unconjugated
Applications	IF,
Format	Liquid
Size	1 ml
Buffer	In 10 mM PBS, pH 7.2 (0.1% sodium azide)
Preservative	0.1% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction S. agalactiae, or GBS, causes pneumonia and meningitis in neonates and the elderly, with

occasional systemic bacteremia. They can also colonize the intestines and the female reproductive tract, increasing the risk for premature rupture of membranes during

Keywords	GBS; Streptococcus Group B
-----------------	----------------------------
