



Anti-S. haemolyticus Polyclonal antibody (DPAB-DC4768)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody reacts only with intact pure streptokinase and with no other proteins so far tested. It reacts with plasminogen-streptokinase complexes.
Target	S. haemolyticus
Immunogen	Native purified streptokinase from Streptococcus haemolyticus.
Source/Host	Sheep
Species Reactivity	S. haemolyticus
Conjugate	Unconjugated
Applications	WB, ELISA, Inhibition Assay,
Format	Liquid
Size	1 ml
Buffer	In serum
Preservative	None
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords S. haemolyticus; Streptococcus haemolyticus streptokinase