



Anti-S. pneumoniae Polyclonal antibody (DPAB-DC4770)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Reacts against all Strep pneumoniae serotypes tested.
Target	S. pneumoniae
Immunogen	Native purified Streptococcus pneumoniae species specific membrane proteins.
Source/Host	Pig
Species Reactivity	S. pneumoniae
Conjugate	Unconjugated
Applications	ELISA, Agg,
Format	Liquid
Size	100 µl
Buffer	In antiserum
Preservative	None
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Streptococcus pneumoniae, or pneumococcus, is a Gram-positive, alpha-hemolytic, aerotolerant, aerobic member of the genus Streptococcus. A significant human pathogenic bacterium, S. pneumoniae was recognized as a major cause of pneumonia in the late 19th
--------------	--

Keywords

S. pneumoniae; Streptococcus pneumoniae
