



Anti-Staphylococcus Protein G Polyclonal antibody (DPAB-DC4773)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Staphylococcus Protein G
Immunogen	Native purified Streptococcus species Protein G.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Staphylococcus
Conjugate	Unconjugated
Applications	sELISA,
Format	Lyophilized
Size	2 ml
Buffer	Lyophilized from 20 mM potassium phosphate buffer, 150 mM NaCl, pH 7.2 (0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C on dry atmosphere. After reconstitution with 2 mL of deionized water, store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction Streptococcus is a genus of spherical Gram-positive bacteria belonging to the phylum Firmicutes and the lactic acid bacteria group. Cellular division occurs along a single axis in

these bacteria, and thus they grow in chains or pairs, hence the name—from

Keywords

Streptococcus Protein G
