



Anti-Staphylococcus SLO Polyclonal antibody (DPAB-DC4774)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody react with Streptolysin O (Streptolysin O of Group A, C, and G origins).
Target	Staphylococcus SLO
Immunogen	Recombinant protein corresponding to hemolytic Streptococcus SLO.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Staphylococcus
Conjugate	Unconjugated
Applications	WB, Inhibition Assay,
Format	Liquid
Size	50 µl
Buffer	In antiserum (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Streptococcus pyogenes is a spherical, Gram-positive bacterium that is the cause of group A
---------------------	---

streptococcal infections. *S. pyogenes* displays streptococcal group A antigen on its cell wall. *S. pyogenes* typically produces large zones of beta-hemolysis (the co

Keywords

Streptococcus SL; Streptococcus SLO
