



Anti-Substance P 4-11 (aa 4-11) polyclonal antibody (DPAB-DC4778)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide (conjugated with thyroglobulin) corresponding to amino acids 4-11 of Substance P 4-11.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rat
Conjugate	Unconjugated
Applications	IHC-P, IHC-Fr,
Format	Lyophilized
Size	50 µl
Buffer	Lyophilized from serum
Preservative	None
Storage	Store at 4°C on dry atmosphere. After reconstitution with 50-100 uL distilled water, store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords Substance P 4-11

GENE INFORMATION

Entrez Gene ID [6863](#)

UniProt ID [P20366](#)
