



Anti-Substance P polyclonal antibody (DPAB-DC4780)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Substance P conjugated with BSA.
Source/Host	Pig
Species Reactivity	Rat
Conjugate	Unconjugated
Applications	IHC-Fr,
Format	Lyophilized
Size	50 µl
Buffer	Lyophilized from serum
Preservative	None
Storage	Store at 4°C on dry atmosphere. After reconstitution with 50-100 uL distilled water, store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	In the field of neuroscience, substance P (SP) is a neuropeptide - a substance that functions as a neurotransmitter and as a neuromodulator. To be specific, substance P is an undecapeptide - a peptide composed of a chain of 11 amino acid residues. It belongs
Keywords	SP; Substance P

GENE INFORMATION

Entrez Gene ID [6863](#)

UniProt ID [P20366](#)
