



# Anti-Sulfabenzamide polyclonal antibody (DPAB-DC4783)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Sulfabenzamide(p) conjugated with BTG.
Source/Host	Sheep
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	EIA,
Format	Liquid
Size	100 µl
Buffer	In 20 mM phosphate buffer, 150 mM NaCl, pH 7.2 (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Introduction	Sulfabenzamide is an antibacterial/antimicrobial. Often used in conjunction with sulfathiazole and sulfacetamide (trade name - Sultrin) as a topical, intravaginal antibacterial preparation.
Keywords	Sulfabenzamide; Sulfabenzamide(p)